

Aquaporin 2 Blocking Peptide peptide (Phospho) (STJ504528)

STJ504528

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-Aquaporin 2 Blocking Peptide is synthetically produced from the 200-240 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25386
Gene Symbol	Aqp2
Uniprot ID	AQP2_RAT
Immunogen	Synthetic peptide taken within amino acid region 200-240 on rat Aquaporin 2 protein. Phosphorylated Serine 218 residue.
Immunogen Region	200-240
Specificity	
Immunogen Sequence	