

**BTK Blocking Peptide peptide (Phospho) (STJ504544)**

STJ504544

**GENERAL INFORMATION**

<b>Product Type</b>	Biomolecules
<b>Short Description</b>	Phospho-BTK Blocking Peptide is synthetically produced from the 544-561 sequence and is suitable for use in western blot applications.
<b>Applications</b>	Immunodepletion/Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

**PRODUCT PROPERTIES**

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	695
<b>Gene Symbol</b>	BTK
<b>Uniprot ID</b>	BTK_HUMAN
<b>Immunogen</b>	Synthetic peptide corresponding to Human BTK amino acids 544-561 (phospho Y551). Sequence RYV LDD EY (p) T SSV GSK FPV
<b>Immunogen Region</b>	544-561
<b>Specificity</b>	
<b>Immunogen Sequence</b>	