

ADCY7 Positive Control for STJ502074 peptide (STJ504558)

STJ504558

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	ADCY7 Positive Control for STJ502074 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	11513
Gene Symbol	Adcy7
Uniprot ID	ADCY7_MOUSE
Immunogen	
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with ADCY7 antibody STJ502074.
Immunogen Sequence	