

AKT1 Positive Control peptide (Phospho) (STJ504563)

STJ504563

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Biomolecules |
| Short Description | Phospho-AKT1 Positive Control is synthetically produced from the sequence and is suitable for use in western blot applications. |
| Applications | WB |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | WB: 1:500 |
| Formulation | |
| Isotype | |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|-------------------------|------------|
| Gene ID | 207 |
| Gene Symbol | AKT1 |
| Uniprot ID | AKT1_HUMAN |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen | |
| Sequence | |