

PAR1 Positive Control for STJ502112 peptide (STJ504569)

STJ504569

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	PAR1 Positive Control for STJ502112 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:250
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	14062
Gene Symbol	F2r
Uniprot ID	PAR1_MOUSE
Immunogen	
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with PAR1 antibody STJ502112.
Immunogen Sequence	