

CACNB2 Blocking Peptide peptide (Phospho) (STJ504581)

STJ504581

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-CACNB2 Blocking Peptide is synthetically produced from the 600-655 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	116600
Gene Symbol	Cacnb2
Uniprot ID	CACB2_RAT
Immunogen	Synthetic peptide taken within amino acid region 600-655 on rat CACNB2 protein.
Immunogen Region	600-655
Specificity	
Immunogen	
Sequence	