

PDE1A Blocking Peptide for STJ502226 peptide (STJ504600)

STJ504600

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	PDE1A Blocking Peptide for STJ502226 is synthetically produced from the 1-50 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide corresponding to unique amino acid sequence within amino acid region 1-50 on PDE1A protein. Common to all PDE1A variants.
Immunogen Region	1-50
Specificity	This blocking peptide is recommended for use in combination with PDE1A antibody, STJ502226
Immunogen Sequence	