

PDE1C Blocking Peptide for STJ502241 peptide (STJ504609)

STJ504609

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	PDE1C Blocking Peptide for STJ502241 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5137
Gene Symbol	PDE1C
Uniprot ID	PDE1C_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on PDE1C gene.
Immunogen Region	
Specificity	This blocking peptide is recommended for use in combination with PDE1C antibody, STJ502241
Immunogen Sequence	