

PDE3A Blocking Peptide for STJ502266 peptide (STJ504612)

STJ504612

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Biomolecules |
| Short Description | PDE3A Blocking Peptide for STJ502266 is synthetically produced from the sequence and is suitable for use in western blot applications. |
| Applications | |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|-----------------------|--|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | |
| Formulation | Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice. |
| Isotype | |
| Storage | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |
| Instruction | |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 5139 |
| Gene Symbol | PDE3A |
| Uniprot ID | PDE3A_HUMAN |
| Immunogen | CQATGYPLAEVSWPNVSPA |
| Immunogen Region | |
| Specificity | This blocking peptide is recommended for use in combination with PDE3A antibody, STJ502266 |
| Immunogen Sequence | |