

PDE3A Blocking Peptide for STJ502266 peptide (STJ504614)

STJ504614

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	PDE3A Blocking Peptide for STJ502266 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5139
Gene Symbol	PDE3A
Uniprot ID	PDE3A_HUMAN
Immunogen	Synthetic peptide selective only for PDE3A variants
Immunogen Region	
Specificity	This blocking peptide is recommended for use in combination with PDE3A antibody, STJ502266
Immunogen Sequence	