

PDE4B1 Blocking Peptide for STJ502306 peptide (STJ504630)

STJ504630

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	PDE4B1 Blocking Peptide for STJ502306 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	24626
Gene Symbol	Pde4b
Uniprot ID	PDE4B_RAT
Immunogen	Synthetic peptide selective only for PDE4B1 subtype variant.
Immunogen Region	
Specificity	This blocking peptide is recommended for use in combination with PDE4B1 antibody, STJ502306
Immunogen Sequence	