

PDE5A Positive Control for STJ502350 peptide (STJ504658)

STJ504658

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	PDE5A Positive Control for STJ502350 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500-1:1, 000
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	8654
Gene Symbol	PDE5A
Uniprot ID	PDE5A_HUMAN
Immunogen	
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with PDE5A antibody STJ502350.
Immunogen Sequence	