

## PDE7A Positive Control for STJ502369 peptide (STJ504667)

STJ504667

### GENERAL INFORMATION

<b>Product Type</b>	Biomolecules
<b>Short Description</b>	PDE7A Positive Control for STJ502369 is synthetically produced from the sequence and is suitable for use in western blot applications.
<b>Applications</b>	WB
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	WB: 1:500
<b>Formulation</b>	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	18583
<b>Gene Symbol</b>	Pde7a
<b>Uniprot ID</b>	PDE7A_MOUSE
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	This is positive control is recommended for use in combination with PDE7A antibody STJ502369.
<b>Immunogen Sequence</b>	