

PDE8A Blocking Peptide for STJ502380 peptide (STJ504675)

STJ504675

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	PDE8A Blocking Peptide for STJ502380 is synthetically produced from the C-Term sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: IP: 1:250 IHC: 1:250 ICC: 1:250 IF: 1:250 CFM: 1:250
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic cyclic peptide from C terminus common to all PDE8A variants.
Immunogen Region	C-Term
Specificity	This blocking peptide is recommended for use in combination with PDE8A antibody, STJ502380
Immunogen Sequence	