

Particulate Guanylate Cyclase C Positive Control for STJ502461 peptide (STJ504703)

STJ504703

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Particulate Guanylate Cyclase C Positive Control for STJ502461 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25711
Gene Symbol	Gucy2c
Uniprot ID	GUC2C_RAT
Immunogen	
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with Particulate Guanylate Cyclase C antibody STJ502461.
Immunogen Sequence	