

## Phantom Blocking Peptide for STJ502490 peptide (STJ504719)

STJ504719

### GENERAL INFORMATION

<b>Product Type</b>	Biomolecules
<b>Short Description</b>	Phantom Blocking Peptide for STJ502490 is synthetically produced from the N-Term sequence and is suitable for use in western blot applications.
<b>Applications</b>	Immunodepletion/Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
<b>Isotype</b>	
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	N-epitope synthetic peptide corresponding to unique amino acid sequences on D. melanogaster Phantom protein.
<b>Immunogen Region</b>	N-Term
<b>Specificity</b>	This blocking peptide is recommended for use in combination with Phantom antibody, STJ502490
<b>Immunogen Sequence</b>	