

Pin1 Blocking Peptide for STJ502496 peptide (STJ504725)

STJ504725

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Pin1 Blocking Peptide for STJ502496 is synthetically produced from the 118-168 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 118-168 on human peptidyl-prolyl cis-trans isomerase NIMA-interacting 1 protein.
Immunogen Region	118-168
Specificity	This blocking peptide is recommended for use in combination with Pin1 antibody, STJ502496
Immunogen Sequence	