

p27 KIP1 Blocking Peptide peptide (Phospho) (STJ504739)

STJ504739

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-p27 KIP1 Blocking Peptide is synthetically produced from the 148-198 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 148-198 on human Cyclin-dependent kinase inhibitor 1 (p27, KIP1) protein.
Immunogen Region	148-198
Specificity	
Immunogen	
Sequence	