

Phospholipase A2 Blocking Peptide peptide (Phospho) (STJ504742)

STJ504742

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-Phospholipase A2 Blocking Peptide is synthetically produced from the 30-80 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9373
Gene Symbol	PLAA
Uniprot ID	PLAP_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 30-80 on human Phospholipase A2-activating protein.
Immunogen	30-80
Region	
Specificity	
Immunogen	
Sequence	