

Phospholipase A2 Positive Control peptide (Phospho) (STJ504743)

STJ504743

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-Phospholipase A2 Positive Control is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9373
Gene Symbol	PLAA
Uniprot ID	PLAP_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	