

MAP3K12 Positive Control peptide (Phospho) (STJ504750)

STJ504750

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-MAP3K12 Positive Control is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7786
Gene Symbol	MAP3K12
Uniprot ID	M3K12_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	