

MERTK Blocking Peptide peptide (Phospho) (STJ504753)

STJ504753

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-MERTK Blocking Peptide is synthetically produced from the 746-757 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	65037
Gene Symbol	Mertk
Uniprot ID	MERTK_RAT
Immunogen	Synthetic peptide corresponding to amino acids 746-757 of Human MERTK. Sequence: CKKIYpSGDYpYpRQG
Immunogen Region	746-757
Specificity	
Immunogen	
Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081