

PMS2 Blocking Peptide for STJ502561 peptide (STJ504755)

STJ504755

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	PMS2 Blocking Peptide for STJ502561 is synthetically produced from the 220-270 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5395
Gene Symbol	PMS2
Uniprot ID	PMS2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 220-270 on human Mismatch Repair endonuclease PMS2 protein.
Immunogen Region	220-270
Specificity	This blocking peptide is recommended for use in combination with PMS2 antibody, STJ502561
Immunogen Sequence	