

PDE3B Blocking Peptide peptide (Phospho) (STJ504779)

STJ504779

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-PDE3B Blocking Peptide is synthetically produced from the 260-310 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5140
Gene Symbol	PDE3B
Uniprot ID	PDE3B_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 260-310 on human PDE3B protein. Selective only for PDE3B variant.
Immunogen Region	260-310
Specificity	
Immunogen	
Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081