

PKC eta Positive Control peptide (Phospho) (STJ504791)

STJ504791

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-PKC eta Positive Control is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5583
Gene Symbol	PRKCH
Uniprot ID	KPCL_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	