

PKC eta Blocking Peptide peptide (Phospho) (STJ504792)

STJ504792

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-PKC eta Blocking Peptide is synthetically produced from the 500-550 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5583
Gene Symbol	PRKCH
Uniprot ID	KPCL_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 500-550 on human PKC eta protein. Phosphorylated T513.
Immunogen Region	500-550
Specificity	
Immunogen	
Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081