

PRP19 Blocking Peptide peptide (Phospho) (STJ504793)

STJ504793

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-PRP19 Blocking Peptide is synthetically produced from the 200-300 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	246216
Gene Symbol	Prpf19
Uniprot ID	PRP19_RAT
Immunogen	Phosphorylated synthetic peptide taken within amino acid region 200-300 on rat PRP19F protein.
Immunogen Region	200-300
Specificity	
Immunogen	
Sequence	