

Prolactin Blocking Peptide peptide (STJ504801)

STJ504801

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Prolactin Blocking Peptide is synthetically produced from the 1-100 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5617
Gene Symbol	PRL
Uniprot ID	PRL_HUMAN
Immunogen	N-epitope specific synthetic peptide taken within amino acid region 1-100 on human Prolactin protein.
Immunogen Region	1-100
Specificity	
Immunogen Sequence	