

Slingshot homolog 1 Blocking Peptide peptide (Phospho) (STJ504816)

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-Slingshot homolog 1 Blocking Peptide is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54434
Gene Symbol	SSH1
Uniprot ID	SSH1_HUMAN
Immunogen	Phosphorylated synthetic peptide corresponding to unique amino acid sequence on Slingshot homolog 1 protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	