

Src Blocking Peptide peptide (Phospho) (STJ504820)

STJ504820

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-Src Blocking Peptide is synthetically produced from the 200-300 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6714
Gene Symbol	SRC
Uniprot ID	SRC_HUMAN
Immunogen	N-epitope phosphorylated synthetic peptide taken within amino acid region 200-300 on human Src protein. Phosphorylated Tyr-215 residue.
Immunogen Region	200-300
Specificity	
Immunogen Sequence	