

STAT3 Blocking Peptide peptide (Phospho) (STJ504823)

STJ504823

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-STAT3 Blocking Peptide is synthetically produced from the 700-800 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6774
Gene Symbol	STAT3
Uniprot ID	STAT3_HUMAN
Immunogen	Phosphorylated synthetic peptide taken within amino acid region 700-800 on human STAT3 protein.
Immunogen Region	700-800
Specificity	
Immunogen	
Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081