

## TORC2 Blocking Peptide peptide (Phospho) (STJ504829)

### GENERAL INFORMATION

<b>Product Type</b>	Biomolecules
<b>Short Description</b>	Phospho-TORC2 Blocking Peptide is synthetically produced from the sequence and is suitable for use in western blot applications.
<b>Applications</b>	Immunodepletion/Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	200186
<b>Gene Symbol</b>	CRTC2
<b>Uniprot ID</b>	CRTC2_HUMAN
<b>Immunogen</b>	Phosphorylated synthetic peptide corresponding to unique amino acid sequence on TORC2 protein.
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	