

## TICAM2 Positive Control peptide (Phospho) (STJ504834)

STJ504834

### GENERAL INFORMATION

<b>Product Type</b>	Biomolecules
<b>Short Description</b>	Phospho-TICAM2 Positive Control is synthetically produced from the sequence and is suitable for use in western blot applications.
<b>Applications</b>	WB
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>
<b>Clone ID</b>
<b>Concentration</b>
<b>Conjugation</b>
<b>Purification</b>
<b>Dilution Range</b>
<b>Formulation</b>
<b>Isotype</b>
<b>Storage Instruction</b>

### TARGET INFORMATION

<b>Gene ID</b>	100302736/353376
<b>Gene Symbol</b>	TICAM2
<b>Uniprot ID</b>	TCAM2_HUMAN
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	