

trp Blocking Peptide peptide (Phospho) (STJ504835)

STJ504835

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-trp Blocking Peptide is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Phosphorylated synthetic peptide corresponding to unique amino acid sequences on D. melanogaster Transient Receptor Potential isoform A protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	