

RCE1 Blocking Peptide for STJ502759 peptide (STJ504848)

STJ504848

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	RCE1 Blocking Peptide for STJ502759 is synthetically produced from the 300-350 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9986
Gene Symbol	RCE1
Uniprot ID	FACE2_HUMAN
Immunogen	Synthetic peptide taken within amino acid sequence 300-350 on human RCE1 protein.
Immunogen Region	300-350
Specificity	This blocking peptide is recommended for use in combination with RCE1 antibody, STJ502759
Immunogen Sequence	