

Robo3 Blocking Peptide for STJ502785 peptide (STJ504856)

STJ504856

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Robo3 Blocking Peptide for STJ502785 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	29546
Gene Symbol	Homer1
Uniprot ID	HOME1_RAT
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on human RIG1 protein.
Immunogen Region	
Specificity	This blocking peptide is recommended for use in combination with Robo3 antibody, STJ502785
Immunogen Sequence	