

Sall4C Blocking Peptide for STJ502872 peptide (STJ504884)

STJ504884

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Sall4C Blocking Peptide for STJ502872 is synthetically produced from the 1-50 sequence and is suitable for use in western blot applications.
Applications	ELISA/IHC/IP/WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human SALL4C protein.
Immunogen Region	1-50
Specificity	This blocking peptide is recommended for use in combination with Sall4C antibody, STJ502872
Immunogen Sequence	