

SGLT1 Blocking Peptide peptide (STJ504919)

STJ504919

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	SGLT1 Blocking Peptide is synthetically produced from the 402-422 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6523
Gene Symbol	SLC5A1
Uniprot ID	SC5A1_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids 402-422 of Human SGLT1.
Immunogen Region	402-422
Specificity	
Immunogen Sequence	