

SIRT1 Blocking Peptide for STJ503029 peptide (STJ504932)

STJ504932

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	SIRT1 Blocking Peptide for STJ503029 is synthetically produced from the 25-75 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	93759
Gene Symbol	Sirt1
Uniprot ID	SIRT1_MOUSE
Immunogen	Synthetic peptide taken within amino acid 25-75 on rat Sirtuin 1 protein.
Immunogen Region	25-75
Specificity	This blocking peptide is recommended for use in combination with SIRT1 antibody, STJ503029
Immunogen Sequence	