

Stra6 Blocking Peptide for STJ503149 peptide (STJ504966)

STJ504966

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Stra6 Blocking Peptide for STJ503149 is synthetically produced from the 200-300 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	20897
Gene Symbol	Stra6
Uniprot ID	STRA6_MOUSE
Immunogen	Synthetic peptide taken within third large intra cytoplasmic loop within amino acid region 200-300 on rat Stra6 protein.
Immunogen Region	200-300
Specificity	This blocking peptide is recommended for use in combination with Stra6 antibody, STJ503149
Immunogen Sequence	