

Syndecan-1 Blocking Peptide peptide (STJ504977)

STJ504977

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Syndecan-1 Blocking Peptide is synthetically produced from the 250-350 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6382
Gene Symbol	SDC1
Uniprot ID	SDC1_HUMAN
Immunogen	C-epitope synthetic peptide taken within amino acid region 250-350 on human Syndecan-1 protein.
Immunogen Region	250-350
Specificity	
Immunogen Sequence	