

TLR3 Blocking Peptide for STJ503277 peptide (STJ505004)

STJ505004

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	TLR3 Blocking Peptide for STJ503277 is synthetically produced from the 650-700 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7098
Gene Symbol	TLR3
Uniprot ID	TLR3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 650-700 on human TLR3 protein.
Immunogen Region	650-700
Specificity	This blocking peptide is recommended for use in combination with TLR3 antibody, STJ503277
Immunogen Sequence	