

TPCN1 Blocking Peptide for STJ503363 peptide (STJ505025)

STJ505025

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	TPCN1 Blocking Peptide for STJ503363 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	53373
Gene Symbol	TPCN1
Uniprot ID	TPC1_HUMAN
Immunogen	Synthetic peptide taken within TMD7 and TMD8 on human TPC1 protein.
Immunogen Region	
Specificity	This blocking peptide is recommended for use in combination with TPCN1 antibody, STJ503363
Immunogen Sequence	