

Visfatin Blocking Peptide peptide (STJ505066)

STJ505066

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Visfatin Blocking Peptide is synthetically produced from the 400-490 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: IHC: 1:100
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	10135
Gene Symbol	NAMPT
Uniprot ID	NAMPT_HUMAN
Immunogen	C-epitope synthetic peptide taken within amino acid region 400-490 on human Visfatin protein
Region	400-490
Specificity	
Immunogen	
Sequence	