

ROCK1 Positive Control for STJ502819 peptide (STJ505102)

STJ505102

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	ROCK1 Positive Control for STJ502819 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution	WB: 1:500
Range	IHC: 1:100
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6093
Gene Symbol	ROCK1
Uniprot ID	ROCK1_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	This is positive control is recommended for use in combination with ROCK1 antibody STJ502819.
Immunogen	
Sequence	