

Soluble Guanylyl Cyclase Alpha 2 Positive Control for STJ502974 peptide (STJ505138)

STJ505138

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Soluble Guanylyl Cyclase Alpha 2 Positive Control for STJ502974 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:750
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	66012
Gene Symbol	Gucy1a2
Uniprot ID	GCYA2_RAT
Immunogen	
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with Soluble Guanylyl Cyclase Alpha 2 antibody STJ502974.
Immunogen Sequence	