

Sodium/Hydrogen Exchanger 3 Antibody Positive Control Positive Control peptide (STJ505160)

STJ505160

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Sodium/Hydrogen Exchanger 3 Antibody Positive Control Positive Control is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemelli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6550
Gene Symbol	SLC9A3
Uniprot ID	SL9A3_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	