

TM6SF2 Positive Control for STJ503312 peptide (STJ505215)

STJ505215

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	TM6SF2 Positive Control for STJ503312 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	107770
Gene Symbol	Tm6sf2
Uniprot ID	TM6S2_MOUSE
Immunogen	
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with TM6SF2 antibody STJ503312.
Immunogen Sequence	