

## TrkA Positive Control peptide (STJ505233)

STJ505233

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Biomolecules  |
| <b>Short Description</b> | TrkA Positive Control is synthetically produced from the sequence and is suitable for use in western blot applications. |
| <b>Applications</b>      | WB  |
| <b>Host/Source</b>       |   |
| <b>Reactivity</b>        |   |

### PRODUCT PROPERTIES

|                       |  |
|-----------------------|--|
| <b>Clonality</b>      |  |
| <b>Clone ID</b>       |  |
| <b>Concentration</b>  |  |
| <b>Conjugation</b>    |  |
| <b>Purification</b>   |  |
| <b>Dilution Range</b> | WB: 1:500  |
| <b>Formulation</b>    | Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is |
| <b>Isotype</b>        |  |
| <b>Storage</b>        | Store at -20°C for long term storage. Avoid freeze-thaw cycles.  |
| <b>Instruction</b>    |  |

### TARGET INFORMATION

|                    |             |
|--------------------|-------------|
| <b>Gene ID</b>     | 4914        |
| <b>Gene Symbol</b> | NTRK1       |
| <b>Uniprot ID</b>  | NTRK1_HUMAN |
| <b>Immunogen</b>   |             |
| <b>Immunogen</b>   |             |
| <b>Region</b>      |             |
| <b>Specificity</b> |             |
| <b>Immunogen</b>   |             |
| <b>Sequence</b>    |             |