

TRPM8 Positive Control for STJ503412 peptide (STJ505235)

STJ505235

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	TRPM8 Positive Control for STJ503412 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	171384
Gene Symbol	Trpm8
Uniprot ID	TRPM8_RAT
Immunogen	
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with TRPM8 antibody STJ503412.
Immunogen Sequence	